



INTEGRATED SECURITY FOR Seamless Care

Every baby needs to be safe, secure, and well cared for. How can hospitals take their security and care to the next level? By using advanced technology, technology gives nurses the flexibility they need to provide seamless care to patients while keeping your facilities protected zones intact—ensuring that safety is never compromised. With caring nurses, skilled doctors, attentive staff, and advanced security technology, your patients are in the best of hands.

- ✓ 4 seconds is all it takes for an abductor to leave a room with a baby and make a hasty exit
- ✓ The majority of U.S. infant abductions occur in hospitals

Seamless care and uninterrupted security

There are many nuances and details when it comes to providing the proper care for patients. But health and safety will always be of utmost importance. Unfortunately, incidents can still occur and infant abductions are still a very real threat. The right balance of security and accessibility will give your patients both the freedom and safety that they need.

What is available for your facility?

You want your facility to be a safe and positive environment for patients. With technology, you get peace-of-mind, safety, and flexibility. What technology is new to the market?

- Stand alone systems that won't go down even if your network does
- Automatic tag transfers that allow nurses to focus more on patient care and less on administrative duties
- Waterproof tags that will allow mothers to shower after giving birth without fear of condensation issues
- Bands that can be used on babies of any size (premes on up)—without false alarms or discomfort to the patient
- Handheld tag testers that allow you to access pertinent tag information on the go





THE TAG identifies each patient as they approach a transfer point (whether it's a door, elevator, or stairway). Make sure your tag is reliable, small, easy to clean, and doesn't come in contact with the baby's skin.

THE BAND holds the tag in place. Make sure your bands are durable, but also comfortable for patients or any size. Bands should also be skin sensing, self-adjusting, and don't cut the baby's skin or hinder circulation.

DOOR CONTROLLERS will respond to the tag on each baby's band. If the patient is cleared to exit, the doors will open. Otherwise an alarm will sound or the door will be locked.

TRANSFER POINTS can be configured based on each patient and can either permit a patient through or not. You can choose to have an alarm sound and lock down the perimeter, different units, or the entire building.

TAG TEST STATION can be used to view remaining battery life, tag notes, and run reports. These stations should be compatible with your facilities current software to enhance efficiency.



Give your patients more freedom and allow nurses more flexibility when transporting patients—all while keeping everyone safe and comfortable.



To learn more about how your facilities could be better protected by software-enhanced Cuddles Infant Protection, contact us today. As an industry leader, we provide high quality, cost-effective, security systems backed by the industry's most comprehensive customer support program with **24/7/365 FREE** technical support and local service.

