


Features you'll only find with the ResidentGuard systems



STAD

The Secure Tag Activator Deactivator (STAD) sets a new standard for cost-effective wander management by lengthening tag life and saving you money. Tags can be tested for low battery using the portable handheld device and unused tags can be turned off using a secure code ensuring long battery life and minimal maintenance.

418 MHz Frequency



All ResidentGuard systems are shielded against RF interference for the life of the system. Each tag operates at a unique frequency (418 MHz) on the edge of a range available only to the Department of Defense. The embedded ID code on the RFID chip in addition to the unusual frequency it is transmitted on combine to make ResidentGuard the most reliable wander management system available—minimizing the problem of false alarms.

Tag



Our long-life tags are securely attached to comfortable, easy-to-clean, water-resistant adjustable wrist bands worn by residents similar to a watch. A red, pulsating LED serves as a continuous visual indication for your staff that the tag is active. The tag can be turned on or off at any time by the STAD to conserve battery power.